



## Detachable and adjustable wheelchair mount

### **BJ-605**

Firm hold for tablets, communicators and other devices

Arm that allows different **equipment like laptops**, electronic planners, joysticks, **communicators**, etc., to be **firmly attached** to the **wheelchair**.

Can be mounted and removed or oriented with 90° turns without the need to adjust any screws, and always returns to the exact same position.

- Compatible with any wheelchair (tubes up to 45mm) Rigid arm that reduces vibrations. Maximum load 4 kg.

When placing the order, indicate the device that that must be compatible, to provide the arm with the appropriate tray. There are different sizes and types.

